



THE TRISTEL TRACEABILITY APP

RECORD, TRACE, ANALYSE AND REPORT WITH
TRISTEL'S NEW TRACEABILITY TOOL



RECORD

OPERATOR, DISINFECTANT,
INSTRUMENT AND PATIENT
INFORMATION



DISINFECT

YOUR DEVICE



TRACE

EVERY SINGLE
DISINFECTATION EVENT



ANALYSE

AND REPORT
ON RECORDED DATA

The Tristel Traceability App is designed to provide traceability and compliance in a smart and efficient way. It comprises a Smart App to record every Tristel disinfection event and a Web Portal for data analysis and reporting.

When used in tandem with a Tristel product such as the Tristel Trio Wipes System, the App records operator, disinfectant (LOT and Expiry) and instrument data, as well as a patient ID. Short videos guide the operator through each disinfection event to ensure compliance with user instructions. Upon completion of a disinfection event, the App produces a digital disinfection record with a unique validation code.

Each record is automatically uploaded to the Disinfection Log on both the App and Web Portal, where data can be searched, filtered and extracted into a CSV file.



INCREASED EFFICIENCY & ACCURACY

The App speeds up manual traceability processes by digitalising it. It makes traceability more efficient and more accurate: it is used while disinfecting instead of afterwards and any handwritten content is eliminated.



VOICE RECOGNITION

The Traceability App is available on a voice command device (VCD). Operators can move through the disinfection event by saying "OK" instead of tapping OK on the smart phone, reducing the risk of cross-contamination.



SECURITY & DATA PROTECTION

The App is hosted by Tristel on a secure Microsoft Azure Cloud server. Microsoft has taken measures to secure data and has more than 20 Cloud computing related security certificates in place, including ISO 27001 (information security) and 27018 (personal data protection). Microsoft also complies with both international and industry-specific compliance standards and participates in rigorous third-party audits that verify security controls.

In addition, the App has been developed according to the latest General Data Protection Regulation, which came into force on 25th May 2018.



PAPERLESS OPERATION

The App is completely paperless. Disinfection records produced by the App are digital. Synchronisation between the App and Web Portal occurs when a WiFi connection is established.



ANALYSIS & REPORTING

The Web Portal gathers all disinfection records in an easy-to-navigate secure environment. The Web Portal is also used to add instruments to be disinfected, operators and clinical settings. This results in an organised overview of data essential for traceability and compliance.

Developed By:

Tristel[™]




GETTING STARTED IS EASY!

TRISTEL ADMIN SETS UP YOUR ACCOUNT, INCLUDING:

- A. ONE OR MULTIPLE SUPER USERS:**
these are members of hospital staff with access to the Web Portal for data analysis and reporting purposes.
- B. CLINICAL SETTING(S):**
the department(s) in which the App will be used.
- C. DEVICE(S):**
which instrument(s) is/are going to be disinfected.
- D. APP USER(S):**
who will be using the App.

UPON RECEIPT OF YOUR TRISTEL APP EQUIPMENT, SIMPLY ENTER YOUR APP PASSWORD TO GET STARTED!



WHAT'S INCLUDED

- Smart phone (VCD and App installed)
- Access to the Web Portal
- Setup Instructions
- Desktop Charger
- Wireless Charger
- Service (3 years)
- Warranty (1 year)
- Extended warranty



AVAILABLE NOW IN:

United Kingdom: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - W www.tristel.com

Ireland: Brennan & Company, 61 Birch Avenue, Stillorgan Industrial Park, Stillorgan, Co. Dublin, A94 XW68 - T 01 2952501 - E enquiries@brennanco.ie - W www.brennanco.ie

Hong Kong: Tristel Asia Ltd, 21st Floor, 168 Electric Road, Fortress Hill, Hong Kong
T +852 2895 6968 - F +852 2869 4388 - E customerservice@tristel.com.hk

Australia: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand: Tristel New Zealand Ltd, 23 Birch Avenue, Judea, Tauranga, P.O. Box 3110 - T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

For any enquiries, please contact your local Tristel representative or email AppSupport@tristel.com